- DE Gebrauchsanweisung Lichtleitkabel
- EN Instruction for Use Fibre Optic Light Cable
- ES Manual de instrucciones de la cable guía de luz
- FR Manuel d'instruction câble de guidage de la lumière
- IT Manuale d'istruzioni de la cavo di guida della luce

A024

DE Diese Gebrauchsanweisung betrifft:
EN These instructions for use concern:
ES Estas instrucciones de uso se refieren a:
FR Ce mode d'emploi concerne:
IT Queste istruzioni per l'uso riguardano:

REF

E.8000.01 H	E.8002.04 H	E.8005.09 H	E.8006.14 H	E.8230.35 LED	E.8239.48 LED	E.8309.35 LED
E.8001.02 H	E.8003.05 H	E.8006.11 H	E.8006.20 H	E.8230.48 LED	E.8300.35 LED	E.8309.48 LED
E.8002.03 H	E.8003.06 H	E.8006.13 H	E.8008.15 H	E.8239.35 LED	E.8300.48 LED	

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Table of contents

Introduction	11
Explanation of symbols.....	11
Warnings.....	12
Intended Use.....	12
Preparation.....	13
Sterilization.....	14
Application.....	14
Maintenance, inspection.....	16
Storage.....	16
Warranty / Repair.....	16
Disposal.....	16
Warranty, service and manufacturer address	17

Introduction

We thank you for your confidence in our products and ask you to observe our instructions. The correct observance of these instructions will give you the guarantee to obtain the optimal function and to allow the longest possible life of the Gimmi® products.

With the purchase of this instrument, you receive a high-quality product whose proper handling and use is described below. In order to minimize risks and avoidable exposure for patients, users and third parties, please read the instructions for use carefully and keep them in a safe place.

Explanation of symbols

Attention! Check the operating instructions for important safety-related information, such as warnings and precautions.

Observe the instructions for use.
Read the instructions for use before application.

Observe the instructions for use.
Read the instructions for use before application.

Manufacturer of the medical device

CE marking for European approval

Reference number

Serial number of the manufacturer.

Non-sterile- Prepare before first use

The instrument must be disposed of separately from the residual waste/private waste system

Caution glass, fragile

Warnings

- Attention! Check the operating instructions for important safety-related information, such as warnings and precautions.
- The medical devices are supplied non-sterile and must be cleaned, disinfected and sterilized before the first use and each subsequent use.
- Defective products must not be used and must have undergone the entire reprocessing process before being returned. A written report to the manufacturer must also be made.
- The light guide cable used does not conduct current due to appropriate insulation. However, if the cold light source used does not provide adequate grounding protection, the metal plug-in socket of the cold light source could be voltages are present which can lead to electric shock. Therefore, always ensure that medically approved cold light sources are used.
- Please observe additional instructions enclosed with the product!
- Remove all protective sleeves and protective films before first use or preparation. As a result, sharp points may be exposed. Handle these areas with the necessary care!
- The safe combination of the products with each other or of the products with implants must be verified by the user before clinical use
- Avoid improperly throwing or dropping instruments.
- After use of the product in patients with Creutzfeldt-Jakob disease or HIV infection, it must be safely disposed of immediately in accordance with applicable national and international guidelines and laws.
- The instrument may only be serviced, repaired or modified by Gimmi® GmbH or persons authorized by the manufacturer.
- Heat is also transmitted with the light guide. In particular, soiling of the light entry or light exit surfaces can lead to increased heating of the light guide connection. Appropriate precautions must therefore be taken (use of gloves or allow to cool down). It is therefore advisable to use a cold light source equipped with a heat protection filter, in which case heat generation at the light guide outlet (instrument side), which could be hazardous to persons, is largely excluded.

Intended Use

Product description

Light guide cables from Gimmi® GmbH make it possible to illuminate the operating area during endoscopic examinations through natural or surgical body openings. Here, the light guide cable serves to transmit the light from the light source to the endoscope, from where the light is then transported into the interior of the body. This allows internal organs, tissues and other structures to be seen, evaluated and treated. The intended duration of use of the instruments is usually less than 24 hours.

General structure

The light guides contain glass fibers as the light-conducting element, consisting of a core, a sheath and a protective coating. The light guide ends are constructed with a metal/plastic combination. Mechanical protection of the fibers is achieved by an inner metal spiral, a braided sleeve and an outer, thick-walled sheath made of a tear-resistant silicone.

To prevent kinking of the hose behind the end pieces, each light guide is additionally reinforced at this point. The reinforcing grommets are fitted on the inside. This facilitates cleaning (no turning upside down when rubbing) and prevents hidden dirt deposits.

Intended Use

Light guide cables are used to transport light from the cold light source to the endoscope to illuminate the surgical field. It is thus required for examining and treating the patient's anatomical structures. The light guide cables are connected to the optics and the light

source by adapters. This enables the transport of light from the light source to the surgical site. Due to the illumination of the application site generated in this way, it can be used for both diagnostic and therapeutic applications.

Indications

Fiber optic cables are used for all visual applications that meet the intended purpose and where the patient requires endoscopic diagnosis and/or therapy and there are no contraindications.

Contraindications

Light guide cables should not be used when endoscopic or microscopic use is contraindicated.

Intended user

The type of application must be determined in each individual case by the surgeon in collaboration with the internist and the anesthesiologist.

The use of fiber optic cables in surgical and/or diagnostic applications in various surgical disciplines must be carried out by appropriately trained and medically qualified personnel.

Before use, the user should carefully read the instructions for use.

Intended patient

Refer to the "Indications" and "Contraindications" items to identify appropriate patients based on the indications and contraindications. Patients of all ages requiring visual diagnosis and/or therapy for whom the fiber optic cable is indicated and on whom its use is necessary and feasible.

Intended place of application

The application environment for these instruments is within a medical facility. Examples include the operating room, ambulances, CSSD, and physicians' offices. Observe the interaction between technical equipment and the instruments. The use of the products must not be impaired or hindered under any circumstances.

Read the instructions for use before use.

These operating instructions must be kept in an easily accessible place near the instrument.

Preparation

1. Manual cleaning

Cleaning brushes with soft bristles can be used as aids for manual cleaning.

The procedure must be followed as follows:

- Immerse the light guide cable completely in the disinfectant and cleaning agent. Make sure that all light guide cables are completely immersed so that they are fully covered with disinfectant and cleaning agent.
- Clean adapter threads at both ends of the light guide cable with a brush.
- Check the light guide cable for residual soiling and brush off again if necessary.
- After disinfecting and cleaning, the light guide cable should be rinsed until all disinfectant and detergent residues have been completely removed.
- The light guide cable ends can still be dried with cleaning rods after the disinfection and cleaning process.

Only brushes that exclude damage to the fiber optic cable may be used.

2. Immersion disinfection

Immersion disinfection is also possible. We recommend the following for immersion disinfection: 2% Sekusept® aktiv (manufacturer Ecolab). Immersion time 30 minutes.

Sterilization

The special design of the light guide cables prevents dirt deposits from developing in the adapter-hose transition. All light guide cables and screw-on adapters can be sterilized in an autoclave at up to 134° C. We recommend an exposure time of 20 minutes. However, sudden evacuation of the light guide cables after sterilization should be avoided in order to prevent possible air pockets in the tube, inflation or even bursting.

Reprocess again if necessary.

To ensure the sterility of the light guide cable, the adapters must be separated from the light guide cable during sterilization.

Machine washing

To remove heavy contamination, it is also possible to clean the light guide cable in a medical washing machine. However, various cleaning agents can damage the silicone tube of the light guide cable! The application profile of the cleaning agent used must be observed.

Application


- Heat is also transmitted with the light guide. In particular, soiling of the light entry or light exit surfaces can lead to increased heating of the light guide connection. Appropriate precautions must therefore be taken (use of gloves or allow to cool down). It is therefore advisable to use a cold light source equipped with a heat protection filter, in which case heat generation at the light guide outlet (instrument side), which could be hazardous to persons, is largely excluded.
- The safe combination of the products with each other or of the products with implants must be verified by the user before clinical use
- The light guide cable used does not conduct current due to appropriate insulation. However, if the cold light source used does not provide sufficient grounding protection, voltages could be present at the metal plug-in socket of the cold light source, which could lead to electric shock. Therefore, always ensure that medically approved cold light sources are used.

Heat and power

When using fiber optic cables, it must always be taken into account that the fiber optic cable is connected and interacts with other products. For this reason, compatibility with surrounding components must always be checked. In addition, the user must take into account that the cable is connected to medical devices. Therefore, if it is not properly connected to the light source or if it is combined with a light source that is not medically approved, it may generate heat or voltage.

Contamination of the light entry or exit surface can also lead to increased heating of the light guide connection.

Using the wrong cable diameter can also lead to increased heating of the connection. For this reason, the cables must be matched to the light entry diameter of the endoscope used. This can prevent damage (burns) to the light exit of the cable caused by reflection. For this purpose, always follow the specification that the diameter of the light guide cable is smaller than the diameter of the endoscope used.

For example, use a \varnothing 3,5 mm optic with a maximum light guide cable diameter of \varnothing 3,5 mm. If, for example, a \varnothing 3 mm optic is used and a light guide cable with a diameter of \varnothing 4,8 mm, this will result in an excess of light and an associated excess of heat, which can lead to heating of the light guide connection and thus to burns.

Mechanical influences

The light guide cable must always be handled with care. Improper handling or damage will hinder optimum application of the cable and may impair the desired result. The matching of the light guide cable to the surrounding instruments, optics and devices must also be carefully thought through.

The fiber optic cable must not be cut open or cut under any circumstances, as the glass fibers it contains can penetrate the skin if they come into contact with it (risk of injury). In addition, liquid can penetrate the cable through the damaged fiber optic cable and interfere with the light transmission. Similarly, germs can settle in cut openings and hinder the effect of disinfection. Do not try to force the plug into the sockets, as there is a risk of injury. The plugs are manufactured to fit exactly, check the type of light guide or adapter if necessary.

Light

The light from the cold light source is guided through the light guide almost without attenuation. Therefore, never look into the exit surface of the light guide. There is a considerable risk of injury to the eye.

Do not look directly into the exit surface of the light guide.

Adaptation

The adaptation to the various cold light sources and instruments is carried out with screw-on adapters. The thread of the cold light cables and screw-on adapters are matched to each other.

The range of different adapters enables the user to combine different cold light sources, instruments and light guide cables from different manufacturers. The adapters are screwed hand-tight onto the threaded connection of the light guide cables (instrument and device side).

Installation

The light guide cable is first fitted with the appropriate screw-on adapter for connection to a cold light source. This is screwed on hand-tight. The light source side of the light guide cable is equipped with a black grip piece. This unit is inserted into the output of the cold light source until it clicks into place. The instrument end, also equipped with adapter for the time being, is adapted to the corresponding instrument. The cable is now ready for operation in conjunction with the cold light source and instrument. After switching on the cold light source, it transmits the light to the instrument. The light guide cable is installed correctly when it is firmly anchored in the adapters.

Handling instructions

The following handling instructions must be observed:

- Pull the cable out of the cold light source only at the handle!
- Do not expose the cable to strong tensile stress!
- Do not bend the cable at too narrow an angle!
- Do not load the cable with heavy objects, run over it or step on it!
- Do not expose the cable unnecessarily to light intensity during work breaks!
- After use, remove the screw-on adapters and clean both components!

All serious incidents occurring in connection with the product must be reported immediately to the manufacturer and the competent authority of the Member State where the user and/or patient is established without delay. Use the following contact addresses on the next page here.

Maintenance, inspection

After cooling down to room temperature, the fiber optic cables must be visually inspected for protein residues and other impurities. Here, slots, locks, closures, tubes and other areas that are difficult to access must be thoroughly inspected. Light guide cables that are not free of residues must be repeatedly subjected to the entire reprocessing process.

Light guide cables that show stains, are bent (light exit or light entry), are no longer functional or are damaged in any other way must be discarded!

Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from -20°C to 60°C. Protect from sunlight and artificial light.

Warranty / Repair

Our products are made of high-quality materials and are carefully checked before delivery. However, even when used as intended, they are subject to wear and tear to a greater or lesser extent depending on the intensity of use. This wear is due to technical reasons and is unavoidable. Should wear-independent faults nevertheless occur, please contact our customer service department. Defective products may no longer be used and must have undergone the entire remanufacturing process before being returned.

All serious incidents occurring in connection with the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established, without delay. Use the following contact addresses here on the next page.

Disposal

Be aware of sharp edges during disposal to prevent injury.

The transposition of European legislation into national law, requires proper disposal of medical devices and instruments. The symbol below indicates that medical products must be disposed of separately from household and residual waste. Please dispose of medical products according to local and national laws.

Dispose of contaminated medical equipment at an appropriate collection point. Medical waste is classified according to the Dangerous Goods Regulation for UN number "UN 3291" (Medical Waste). Please mark the waste container with the symbol below. Dispose of medical waste according to local and national law.

You can dispose of uncontaminated packaging material according to local and national recycling regulations.

Warranty, service and manufacturer address

Gimmi® GmbH provides a warranty on these products, which is recorded within the Gimmi® GmbH warranty guidelines. The warranty only applies to defects that are not caused by normal wear and tear, misuse, incorrect use, overloading and improper handling. In case of repairs or maintenance, contact the Gimmi® GmbH service or an authorized partner:

Gimmi® GmbH

Carl-Zeiss-Str. 6

78532 Tuttlingen / Germany

Tel.: +49 (0)7461 965 90 - 0

Fax: +49 (0)7461 965 90 - 33

Email: info@gimmi.de

www.gimmi.de

In the interest of expediting your request, please send the product with the following information:

- Part number (REF)
- Serial number (SN)
- Detailed error report

Specialty Retailer:

Manufacturer:

Gimmi® GmbH

Carl-Zeiss-Str. 6

78532 Tuttlingen / Germany

Tel.: +49 (0)7461 965 90-0

Fax.: +49 (0)7461 965 90-33

E-mail: info@gimmi.de

www.gimmi.de

If you have any further questions, please feel free to contact us. We are happy to help!

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