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
Table of contents

Explanation of Symbols	10
Contraindications and use restrictions.....	11
Application and safety notes.....	11
Reprocessing.....	11
Preparation at the point of use.....	12
Transport	12
Preparation for decontamination / manual pre-cleaning	12
Automated cleaning.....	12
Visual inspection and servicing.....	13
Packaging	14
Sterilization.....	14
Storage of sterile products	14
Disposal.....	14
Information on the equipment and materials used.....	15
Studies to prove the preparation parameters	15

We thank you for your confidence in our products and ask you to observe our information. Correct compliance with these instructions will ensure that you will receive optimum performance and the longest possible life for Gimmi® trays.

Explanation of Symbols

Caution! Refer to user manual for important safety-related information such as warnings and precautions.

Manufacturer of the medical device

CE marking for European market access

Intended purpose / indication

The trays are used in the reprocessing of medical products. They are used as containers for medical products which have to be sterilized after cleaning and disinfection. It is the user's responsibility to select the appropriate screen basket for the intended use.

Contraindications and use restrictions

The trays are to be used exclusively for cleaning, disinfection, sterilization and storage of medical products. There are no other special contraindications or restrictions on use.

Application and safety notes

- The application as well as the reprocessing (cleaning, disinfection and sterilization) may only be carried out by trained professionals.
- Please note that our trays are supplied non-sterile and must be prepared (cleaned, disinfected and sterilized) before first use.
- If serious incidents occur in connection with the product, these must be communicated to the manufacturer and the competent authorities in the member state in which the user and/or patient is located.
- Do not use a damaged or defective product. Immediately separate out, label, and exclude damaged products from further use.
- Do not overload the trays during application. Overloading by applying too much force, such as overloading or rough handling of the screen basket during the cleaning process, can lead to damage which can adversely affect further use and create a potential source of injury for the user.

Reprocessing

Frequent reprocessing of trays has no effect on the lifespan of the products. The lifespan is normally determined by wear and damage during use.

Risk of infection

- Completely reprocess the product before the first use.
- Reprocess the product before returning to the manufacturer.
- The instructions in the user manuals of the cleaning agents and disinfectants used as well as that of the washers and sterilizers used must be followed.

- Personal protective equipment must be worn during reprocessing (Gloves, eye protection, etc.).
- Special requirements for reprocessing must be met if there is a suspicion of prions and Creutzfeldt-Jakob disease. These procedures are not described in this reprocessing manual.

Preparation at the point of use

To prevent residues from drying and to remove coarse dirt, the following steps must be performed.

- Rinse trays with cold water.
- Do not use fixing agents or hot water ($>40^{\circ}\text{C}$).
- Remove rough dirt with cold water.
- Rinse all joints and cavities thoroughly with cold water.

Transport

Safe storage and transport, e.g. in a closed container to avoid damage to the trays and contamination from the environment.

Preparation for decontamination / manual pre-cleaning

- Place trays in cold water for at least 5 minutes.
- Clean under cold water with a soft brush (do not use metal brushes), open and closed, until your residues are more visible.
- In the case of cavities, bores and threads, pressure flush for at least 10 seconds with a cleaning gun (pulsed method). We recommend the Gimmi® cleaning gun T.1003.00.

Automated cleaning

Perform automated thermal disinfection taking into account the national requirements regarding A_0 values (see ISO 15883).

Machine clean all trays with alkaline cleaners and thermal disinfection. To optimize the cleaning process, we recommend the use of distilled, demineralized or fully demineralized water.

Cleaning	Dosage	Soaking time
Neodisher FA – Dr. Weigert (alkaline)	0,5 %	According to the manufacturer's instructions

Place the screen basket in the washer-disinfector after preliminary cleaning. Follow the instructions of the device manufacturer.

Start cleaning program:

1. Pre-rinse with cold water: 1 min
2. Empty
3. Repeat pre-rinse with cold water: 3 min
4. Empty
5. Clean with 0.5% alkaline cleaning agent for 5 min at 55 °C (follow the directions of the manufacturer of the cleaning agent)
6. Empty
7. Neutralize with deionized water: 3 min
8. Empty
9. Rinse with deionized water: 2 min
10. Empty
11. Dry at 90–110 °C for 15–25 min. The washer-disinfector program must include an appropriate drying phase
12. Open the machine and allow the remaining steam to escape.
13. Dry hollow spaces of instruments with sterile compressed air.

At the end of the program remove the product from the washer.

Visual inspection and servicing

Reliable sterilization is possible only with clean medical products. For this reason, cleaning efficiency must be ascertained by an intensive visual check of all surfaces. The cleaning process must be repeated if remaining contamination, such as encrustations or deposits, is visible.

Correct servicing of the trays extends their lifespan and hence, should be undertaken after each cleaning. After each cleaning and before sterilization, the moving parts and sliding surfaces must be lubricated and maintained with a medicated oil (Gimmi article numbers: J.8860.o6 and J.8860.o8). The oil used (e.g., paraffins in accordance with Ph. Eur.) must not affect the result of the subsequent sterilization.

Packaging

Standard-compliant packaging of trays for sterilization according to ISO 11607 and EN 868.

Sterilization

Sterilization must be carried out in accordance with ISO 17665, taking into account country-specific national requirements.

Where prions and Creutzfeldt-Jakob disease are suspected, special processing requirements must be observed.

Sterilization is performed with 4 x pre-vacuum phases, 134°C and 5 min. holding time (full cycle) at 3 bar (44 psi) pressure.

The following countries have differing specifications that must be met:

Country	Sterilization time
Germany	≥ 5 – 30 min.
France	≥ 18 – 30 min.
Switzerland	≥ 18 – 30 min.

Place trays in the sterilizer so that no components touch each other and the steam can circulate freely.

The drying time is 10 minutes.

Storage of sterile products

The following storage conditions must be maintained:

- Store the sealed sterile product in a clean, dust-free, and dry environment.
- Protect from direct light.
- Protect from high temperature fluctuations. (Storage temperature between + 0°C and + 50°C)
- Do not store sterile products in the vicinity of aggressive substances (e.g., alcohols, acids, bases, solvents, and disinfectants).

Note: Your in-house storage standards for sterilized products must also be met.

Disposal

Valuable raw materials can be recovered through environmentally responsible disposal. Dispose of the product in an environmentally friendly manner in

accordance with the applicable hospital guidelines and country-specific regulations and laws.

- Microbiological or infection risks caused by products contaminated with potentially infectious substances of human origin.
- Risk of injury caused by sharp tips or edges. Dispose of using suitable protective measures.

Information on the equipment and materials used

The following materials and machines were used in the validation:

Cleaning agent:

Neodisher FA from Dr. Weigert (alkaline cleaner)

Washer-disinfector:

G 7735 CD from Miele

Sterilizer:

Autoklav MMM 6-6-6 Selectomat HP

Sterilizing agent:

Steam, 134°C

The user must ensure that the reprocessing process, including resources, material and personnel, is suitable for achieving the required results. The state-of-art and national laws require compliance with the validated processes. If the chemicals and machines described above are not available, the user must validate his procedure accordingly.

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© Copyright Gimmi® GmbH · Tuttlingen · Germany
Printed in Germany · All rights reserved ·
Gimmi® GmbH